


Installation Qualification, Operational Qualification and Requalification for Table-Top Steam Sterilizers

IPAC Program Template

Infection Prevention and Control (IPAC) Program

Last updated January 2025

© 2025 College of Physicians & Surgeons of Alberta



PURPOSE

This document outlines the process for conducting installation qualification (IQ), operational qualification (OQ) and requalification (RQ) of your medical clinic's table-top steam sterilizer to ensure compliance with [CPSA Reusable & Single-Use Medical Device Requirements for Medical Clinics](#) and safe operation.

For assistance or more information, please contact the CPSA IPAC Program at ipac@cpsa.ab.ca or 780-969-5004.

QUALIFICATION/REQUALIFICATION FREQUENCY

Installation qualification (IQ)

- Must perform IQ immediately after installation or if the sterilizer is moved. Please refer to section 8.3 of [CPSA Reusable & Single-Use Medical Device Requirements for Medical Clinics](#) for IQ testing requirements.

Operational qualification (OQ)

- Must perform OQ after installation qualification (IQ) and prior to use. Please refer to section 8.4. of [CPSA Reusable & Single-Use Medical Device Requirements for Medical Clinics](#) for OQ testing requirements.

Requalification (RQ)

- Must perform and document RQ annually. Please refer to section 8.4 of [CPSA Reusable & Single-Use Medical Device Requirements for Medical Clinics](#) for RQ testing requirements.
- Additional triggers for requalification include:
 - Following any significant repair.
 - Sterilizer relocation.
 - Unexplained sterility failure.
 - Any disruption to steam supply or change to steam pressures.

Installation qualification (IQ)

General information

Clinic name: _____

Sterilizer model and serial number: _____

Location of sterilizer: _____

Date of installation qualification: _____

Installation qualification (IQ) process

- Confirm the sterilizer is installed according to the manufacturer's instructions for use (MIFU).
- Verify all connections (e.g., power, water supply) are secure.
- Verify clearances.
- Verify orientation and appropriate surface.
- Complete installation qualification sign-off.

Documentation and records

- Attach all relevant and supporting documentation to this form. Retain records.

Installation qualification sign-off

Performed by: _____

Name	Signature
------	-----------

Date: _____

References

- Sterilizer Manufacturer's Instructions for Use (MIFU)
- [CPSA Reusable & Single-Use Medical Device Requirements for Medical Clinics](#)
- [CPSA Qualification of Table-Top Sterilizers](#)

Operational qualification (OQ) and requalification (RQ)

General information

Clinic name: _____

Sterilizer model and serial number: _____

Location of sterilizer: _____

Date of operational qualification and requalification: _____

The following tests must be conducted and the results documented:

1. Dynamic air removal test/Bowie Dick (if applicable)

- Purpose: To verify air removal and steam penetration.
- Applicability: Only for pre-vacuum sterilizers.
- Criteria: Test sheet must show uniform colour change with no unprocessed areas.

Results: _____

2. Biological indicator (BI) process challenge device (PCD) test

- Purpose: To verify sterilization efficacy.
- Procedure: Conduct three consecutive cycles using a BI PCD and record results and ensure the chamber is fully loaded during each test.
- Criteria: No growth observed in the biological indicators.

Test 1

Cycle number: _____

Load contents: _____

Test 1 result: _____

Test 2

Cycle number: _____

Load contents: _____

Test 2 result: _____

Test 3

Cycle number: _____

Load contents: _____

Test 3 result: _____

Control BI test result: _____

3. Physical parameters verification

- Purpose: To confirm the sterilizer operates within specified temperature, pressure and time ranges.
- Documentation: Record cycle parameters from sterilizer printouts.

Test 1

Cycle temperature: _____

Cycle time: _____

Test 2

Cycle temperature: _____

Cycle time: _____

Test 3

Cycle temperature: _____

Cycle time: _____

4. Chemical indicator test

- Purpose: To ensure the parameters to achieve sterilization have been met.
- Criteria: All indicators must pass as outlined in the indicator MIFU.

Test 1

Results: _____

Test 2

Results: _____

Test 3

Results: _____

Qualification/requalification acceptance criteria

The sterilizer is considered operationally qualified/requalified if all the following conditions are met:

- Uniform colour change in the Bowie-Dick test (if applicable).
- Three consecutive passing BI PCD test results with fully loaded chambers.
- Physical parameters are within acceptable ranges.
- All chemical indicators demonstrate the required changes.

Corrective actions

If any tests fail, the following actions should be taken, at minimum:

- Identify the cause of failure.
- Perform maintenance or calibration as needed.
- Retest and document results.
- Record corrective actions taken: _____

Documentation and records

- Attach all relevant test results, printouts and logs to this form. Retain records.
-

Qualification/requalification sign-off

Performed by: _____
Name Signature

Date: _____

References

- Sterilizer Manufacturer's Instructions for Use (MIFU)
- [CPSA Reusable & Single-Use Medical Device Requirements for Medical Clinics](#)
- [CPSA Qualification of Table-Top Sterilizers](#)